

ISOL/8[®]

SPECIFICATION SHEET

*Designed by Healthcare Professionals,
For Healthcare Professionals*



Comfort Theatre Cap (elasticated back)

File no.: SPEC-IS8HC03

ISOL8 HEALTHCARE LTD

Unit 4 Joyce House,
Barrack Square, Ballincollig, Co Cork,
P31 HW35,
Ireland


CONTACT

info@isol8health.com
IE: +353 (0) 66 71 21659
GB: +44 (0) 191 816 1966
www.isol8health.com


Lightweight, extra soft and breathable material

Our range of Comfort headwear are all made from a lightweight, extra soft spunlace fabric which offers the highest level of comfort to the wearer. Air vents in the fabric provide excellent breathability, which along with its superior comfort makes these hats perfect for prolonged use in the operating theatre.

This single-use comfort theatre cap has an elasticated rear which is perfect for tucking under tied-back and long hair.



Comfort



Comfort

Dimensions

Crown (Length x Width)	22 x 21 cm \pm 1cm
Depth (front)	13cm \pm 1cm
Length/circumference (elastic stretched)	64cm \pm 3cm

Product Specification

Product Name	Comfort Theatre Cap (elasticated back)
Manufacturer Product Codes (MPN)	IS8HC03
NPC	BWF85020
Manufacturer	ISOL8 Healthcare Ltd.
Product type	Single use theatre cap with elasticated back suitable for long hair – standard depth
Sizes (see dimensions below)	One size fits all
Fabric	Soft Spunlace 25gsm (viscose/polyester)
Colour	Blue with paw-print pattern
Sterility	Non-sterile
Ties	Elasticated at the rear
Latex status	Not made with natural rubber latex
Metal status	Nickel and cadmium free
Shelf-life	3 years
Product packaging (inner box)	100 hats per dispenser box 330x140x185 mm
Product packaging (outer carton)	6 boxes per carton (600 hats)
Recycling Information	Outer carton and inner box are 100% recyclable (cardboard) Hats are non-recyclable
CE certification	MDR 2017/745 – Class I medical device
Manufacturer QMS	ISO 13485
GTIN code inner box	05391541270379
GTIN code outer carton	05391541270942



ISOL8 theatre headwear products are CE class I medical devices intended for use by medical professionals. Before use, please consult the product labels and instructions.

For more information or any enquiries regarding our products, please contact your account manager at ISOL8 Healthcare.

© 2023 ISOL8 HEALTHCARE LTD
Unit 4 Joyce House, Barrack Square, Ballincollig,
Co Cork, P31 HW35. Company No. 687076

info@isol8health.com
IE: +353 (0) 66 71 21659
GB: +44 (0) 191 816 1966
www.isol8health.com

ISOL/8[®]

**Designed by Healthcare Professionals,
For Healthcare Professionals**

Date of issue: 22/03/23, Version 002