

COMFORT RANGE

# SUSTAINABLE HEADWEAR



## 100% renewable

Made from fast-growing eucalyptus and birch-pulp



## 100% biodegradable (compostable)

Product can be disposed of as compostable waste



## Softest fabric on the market

Breathable and extra-soft cotton feel




## Up to 50x less emissions

Compared to a home-washed reusable program<sup>1</sup>

### Background

Our previous comfort fabric, used by many companies was composed of 68% viscose and 32% polymer based binders. These polymer based binders are non-biodegradable and toxic to the environment. We therefore worked with leading textiles institutions to develop a method to eliminate these toxic and non-biodegradable binders to create an extra soft 100% viscose fabric made from eucalyptus and birch pulp. Our new fabric is now 100% biodegradable, 100% renewable and non-toxic to the environment. Eucalyptus and Birch are fast growing trees which are harvested and replanted to ensure an ongoing sustainable supply chain.



Our comfort fabric is designed to be extra soft and breathable, suitable for prolonged use in operating theatres and other hospital settings.

### ORDER INFORMATION

Full range on the NHSSC Catalogue

Old fabric (left) vs. new fabric (right).



### Theatre Cap (tie-back)

Most suitable for persons with shorter hair.

- 100 hats/box
- 6 boxes/carton

REF: IS8HC01 NPC: BWF85017




### Theatre Cap (elasticated back)

Elasticated rear for tucking under tied-back and long hair.

- 100 hats/box
- 6 boxes/carton

REF: IS8HC03 NPC: BWF85020



### Operating Cap (tie-back)

Extended depth for greater coverage for persons with longer hair.

- 100 hats/box
- 6 boxes/carton

REF: IS8HC02 NPC: BWF85015



### Operating Cap (elasticated back)

Elasticated rear and extended depth for greater coverage.

- 100 hats/box
- 6 boxes/carton

REF: IS8HC04 NPC: BWF85016



### Surgeon's hood (ties)

Maximum coverage for long procedures and persons with long hair.

- 100 hats/box
- 6 boxes/carton

REF: IS8HC05 NPC: BWF85014

# ISOL/8<sup>®</sup>

Designed by Healthcare Professionals,  
For Healthcare Professionals

sales@isol8health.com  
GB: +44 (0) 191 816 1966  
www.isol8health.com

[1] Based on washing/drying three caps at 70-90 °C as per infection control guidelines